

## IMPLANTATION OF MAGICTOUCH TECHNOLOGY

### 1. LESION



### 2. PRE-DILATATION



**PRE-DILATATION RECOMMENDED IN ALL CASES**  
For pre-dilatation in all cases use a standard balloon (approx 0.5mm smaller than RVD). Choose a MagicTouch SCB with normal size equal to reference diameter.

### 3. MAGICTOUCH SCB



Duration time 30-60 sec



### DRUG RELEASE WITHIN 60 SECONDS

Longer inflation times are possible at discretion of the operator to pursue optimal mechanical dilatation – but this has no effect on additional drug release.


### SINGLE MAGICTOUCH SCB

Each lesion should be addressed with separate balloon



### OVERLAPPING OF MAGICTOUCH SCB

In longer lesions MagicTouch SCB overlapping is indicated




LESION 1, MAGICTOUCH SCB 1



LESION 2, MAGICTOUCH SCB 2

### IN CASE OF ISR

If stent is present, size of balloon should be greater than the stent by 2 mm on each ends.



## TECHNICAL SPECIFICATIONS

### DRUG / EXCIPIENT

DRUG	Sirolimus
DRUG DOSE	1.27 µg/mm <sup>2</sup>
DRUG CARRIER	Phospholipid Based Excipient

### BALLOON

BALLOON MATERIAL	Polyamide
CATHETER DESIGN	Rapid Exchange (RX) Design

### DELIVERY SYSTEM

Shaft Diameter - Proximal	1.7 F	Rated Burst Pressure	16 bar (14 bar for 4.00 / 25 to 40 mm)
Shaft Diameter - Distal	2.5 F	Guiding Catheter Compatibility	5F
Usable Catheter Length	140 cm	Guidewire Compatibility	0.014" maximum recommended
Tip Profile	0.016"		
Nominal Pressure	6 bar		

### ORDERING INFORMATION

Diameter/Length	10 mm	15 mm	20 mm	25 mm	30 mm	35 mm	40 mm
1.50 mm	CMT15010	CMT15015	CMT15020	CMT15025	CMT15030	CMT15035	CMT15040
2.00 mm	CMT20010	CMT20015	CMT20020	CMT20025	CMT20030	CMT20035	CMT20040
2.25 mm	CMT22510	CMT22515	CMT22520	CMT22525	CMT22530	CMT22535	CMT22540
2.50 mm	CMT25010	CMT25015	CMT25020	CMT25025	CMT25030	CMT25035	CMT25040
2.75 mm	CMT27510	CMT27515	CMT27520	CMT27525	CMT27530	CMT27535	CMT27540
3.00 mm	CMT30010	CMT30015	CMT30020	CMT30025	CMT30030	CMT30035	CMT30040
3.50 mm	CMT35010	CMT35015	CMT35020	CMT35025	CMT35030	CMT35035	CMT35040
4.00 mm	CMT40010	CMT40015	CMT40020	CMT40025	CMT40030	CMT40035	CMT40040

The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and user instructions can be found in the product labelling / IFU supplied with each device. For restricted use only in countries where product is registered with applicable health authorities. All cited trademarks are the property of their respective owners.


Concept  
Medical

Concept Medical Inc. • 5600 Mariner ST STE 200, Tampa, FLORIDA 33609-3417.  
Concept Medical B.V. • Hogebrinkerweg 33, 3871KM Hoevelaken, The Netherlands.  
Concept Medical • 901~903, Trinnity Orion, Beside Jolly Residency, Vesu, Surat-395007, Guj-India.  
www.conceptmedical.com info@conceptmedical.com f t @ in /conceptmedicals

MTVER.01.19.05.01

# MagicTouch

## Sirolimus Coated Balloon




# MagicTouch

MagicTouch is a Sirolimus Drug Coated Balloon Catheter for the treatment of coronary artery disease; intended for In-Stent Restenosis, Small Vessels, Bifurcation lesions and De-Novo lesions.

## NANOLUTE TECHNOLOGY

- Conversion of Sirolimus drug into submicron sized particles.
- Encapsulation of submicron sized Sirolimus drug into highly biocompatible drug carrier – Phospholipid.
- Phospholipid comprises of one hydrophilic head and two lipophilic tails, which improves adhesion property of encapsulated Sirolimus.
- Upon inflation of MagicTouch SCB at target site, drug carrier with Sirolimus drug inside gets transferred to the vessel wall following the principle of co-efficient diffusion.
- Upon body PH variation, drug carrier mimics the body lipids and liberates Sirolimus.
- The submicron sized Sirolimus drug particles penetrate the deepest layer of the vessel over a period.




## ADVANTAGES OF NANOLUTE TECHNOLOGY

- Facilitates better adhesion of Sirolimus on the balloon surface
- Effective drug transfer to the deepest layer of the vessel
- Circumferential coating
- Better in-tissue bioavailability of Sirolimus


## DRUG

Drug: Sirolimus  
 Dose: 1.27µg/mm<sup>2</sup>  
 Nature: Less Lipophilic  
 Action: Cytostatic  
 Therapeutic Range: Wide



PHARMACOKINETIC STUDY

## CARRIER




Carrier: Phospholipid  
 Nature: Biocompatible  
 Action: Excipient + Stabilizer

## BALLOON

Balloon: Semi-compliant Recommended  
 Pre-Dilatation Crossing Profile: 0.029"  
 (Vary with size)  
 Entry Profile: 0.016"

## SIROLIMUS DISTRIBUTION STUDY

DTF labelled Sirolimus was used to study the drug distribution following DCB treatment\*



\*EuroIntervention. 2013 May 20;9(1): 148-56

## UNIQUE COATING TECHNOLOGY

